

# **MTGR integrity test**

## **Multiple positions manifold**

Test simultaneously 5 to 10 filters with the alcohol based bubble point test.

### **Testing approach validity**

Parallel testing of aervent 50 capsules:

Testing multiple aervent 50 on a manifold with the bubble point test is an acceptable practice since it is identical to testing multiple filter cartridges in a multi-round housing. This is a standard industry practice, well documented in the PDA technical reports 26 and 40 on sterilizing filtration of liquids and gas.

The integrity test must ensure that :

- The test is sensitive to detect a damaged filter among the other ones.
- The test is accurate.
- The filter traceability is ensured.

Discussion:

1. The sensitivity is ensured by the bubble point test which is a measurement of the largest pore/smallest defect of the system and is independent to the system size.

(minimum bubble point specification is 1.1 bar with IPA/water 60/40 mixture)

A bubble point measure of 1.5 bar due to one filter reaching the bubble point at this pressure is a guaranty that the filter is integral and that all other filters have a bubble point higher than 1.5 bar and are integral.

A bubble point measure of 0.5 bar due to one defective filter should trigger the testing of the filters individually to identify the failing one.

This sensitivity test is easily verified during the IQOQ of the tester and the manifold.

2. Test accuracy is verified during the on-site qualification of the automatic instrument by the comparison of the test in parallel with a reference method (visual detection + reference pressure gauge).

3. Filter traceability is ensured by organizing the flow of vent capsules in the department, and by their identification with the unique combination of lot number and serial number. documenting different serial number on the same test report is acceptable.

### **With hose barb connection**



**With sanitary TC connection**

$\frac{3}{4}$ " TC Manifold allows for the testing of 50mm air filters with the following configurations:

**$\frac{3}{4}$ " TC**

MTGR15000  
MTGR25000



**Luer**

MTGR85000  
MTGR05000



**Hose**

MTGR85000  
MTGR05000



**$\frac{3}{4}$ " TC**

Stopper




## Utilization

### Filter wetting

#### With a pump or a pressure dispensing vessel


With IPA water mixture 70/30 in parts

Pressure dispensing vessel  
or pump



#### Manually with a syringe


Wet individually the Aervent 50 with a syringe prior to connect to the manifold



## Integrity test


regulated  
gas supply

Integrity  
tester




## Qualification

Simultaneous automatic and manual filter test set-up for qualification purpose



## References

Manifold drawing



Manifold Prince : about 1500 €

## Accessories

Reference	Description
JMSZ50553	as described in DB (3/4" TC ISO DN8)
FTPFO1893	Clamp 3/4" TC - stainless steel
or 3538200WB	Clamp 3/4" TC - plastic
MS0TC3403	Adaptator for sampling port TC 3/4"-10 pieces
to be created	Gaskets EPDM
to be created	3/4" TC blindcap
Total 1	stainless steel
Total 2	plastic clamp