



Associates of Cape Cod, Inc.

Your Endotoxin & Glucan Experts

CELEBRATING OUR 45 YEAR ANNIVERSARY

GLUCAN DETECTION



**Our Test, Your Cure...
Ensuring A Healthy World**

Your Endotoxin Experts!

Experts in Endotoxin and Glucan Detection for 45 Years.



GLUCAN DETECTION



Glucatell® Kit – (1→3)-β-D-Glucan Detection

The Glucatell® kit is specific for detection of (1→3)-β-D-Glucan (BG). BG is considered a contaminant in parenteral drugs, solutions, and medical devices due to BG's bioactivity as an activator of innate immunity and its capacity to enhance bacterial endotoxin tests (literature references to the bioactivity of BG may be found at www.acciusa.com/pdfs/BG_Biol_Activity_Ref_List_PR16023.pdf). The Glucatell assay is based upon a modification of the *Limulus* Amebocyte Lysate (LAL) pathway. Glucatell® reagent is processed to eliminate Factor C, and is therefore specific for (1→3)-β-D-Glucan. The reagent does not react with other polysaccharides, including beta-glucans with different glycosidic linkages. Glucatell® is a chromogenic reagent that may be used to perform either kinetic or endpoint assays in microplate readers.

Sensitivity

The sensitivity for Glucatell® assay is determined by the lowest standard concentration on the standard curve used for the assay. The maximum sensitivity for Glucatell® is 3.125 pg/mL when used with a microplate reader.

Sample to Glucatell® Ratio

Kinetic Assay: Reconstituted Glucatell® reagent is used at a ratio of 1:4 and a volume of 100 µL/well (55 tests/vial).

Endpoint Assay: Reconstituted Glucatell® reagent is used at a ratio of 1:1 and a volume of 50 µL/well (55 tests/vial).

Performing the Test

Kinetic Assay: The Glucatell®/sample mixture is incubated at 37±1°C in a microplate reader. Software is used to construct the standard curve and calculate glucan concentrations.

Endpoint Diazo Assay: The Glucatell®/sample mixture is incubated at 37±1°C in a microplate heating block for the recommended time period. 50 µL each of the three diazo reagents are then added to the mixture. Software is used to construct the standard curve and calculate glucan concentrations.

Reconstitution

Glucatell® reagent can be reconstituted differently depending on the assay you use. For endpoint assays, use 2.8 mL of only Pyrosol® reconstitution buffer. For kinetic assays, combine 2.8 mL each of Pyrosol® and Reagent Grade Water.

Stability

Store all reagents at 2–8°C in the dark. Once reconstituted, Glucatell® reagent should be stored at 2–8°C and used within 2 hours. Alternatively, reconstituted Glucatell® reagent can be frozen at -20°C for 20 days, thawed once and used. The diazo reagents should be used the day they are prepared.

Product Applications

- Analyzing final products for (1→3)-β-D-Glucan
- Investigating BET Out of Specification results
- Qualifying raw materials
- Monitoring cellulosic filter extractables
- Monitoring fungal fermentation processes
- Analyzing fermentation and cell culture media
- Monitoring airborne glucan burden

Packaging

The GlucateLL® kit is available as either an endpoint or kinetic chromogenic assay for use in microplates. The kit contains the GlucateLL® reagent, a (1→3)-β-D-Glucan standard, buffer, glucan-free water, glucan-free microplates and diazo reagents (endpoint kit only).

GlucateLL® Kit

- | | |
|--------|---------------------------------------------------|
| #GT002 | Kinetic assays; 110 tests |
| #GT003 | With diazo reagents for endpoint assays; 55 tests |
| #GT004 | Kinetic assays; 55 tests |

Accessory Products

Pipette tips

- | | |
|--------|-----------------------------------------|
| #PPT25 | 250 µL tips, 96 tips/box, 10 boxes/pack |
| #PPT10 | 1000 µL tips, 96 tips/box, 8 boxes/pack |

Pyrotubes®

- | | |
|----------|-----------------------------------------------------------------------------------------|
| #TB013-5 | 13 x 100 mm borosilicate glass
(for dilutions only) 18 tubes/pack,
5 packs/carton |
|----------|-----------------------------------------------------------------------------------------|



ORDERING INFO AND SALES OFFICES



Ordering Information

Customer service representatives are available to assist you with orders, pricing requests and Certificates of Analysis.

Method of Payment For United States

- Check (in US dollars) made payable to Associates of Cape Cod, Inc.
- Wire Transfer (contact Accounts Receivable for routing information)
- Credit Card (AMEX®, VISA®, MasterCard®)
If payment is to be made by credit card, the following information is required: Type of Credit Card, Card Number, Credit Card Security Code, Expiration Date of Card, and Name (as it appears on the card).

Additional Information

ACC reserves the right to institute, modify or discontinue credit limits provided to customers at any time for any or no reason.

The use of credit cards for payment may incur a fee, please see our website for ACC's policy on credit card usage at www.acciusa.com/pdfs/acctc.pdf

Outside the U.S.

Please contact your local office for information regarding method of payment. For a list of your country specific distributors please visit www.acciusa.com.

Product listings, information and fill sizes are subject to change at any time without prior notice.

Corporate Headquarters

Associates of Cape Cod, Inc.

124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536-4445 USA
Tel: (800) LAL-TEST (525-8378) or (508) 540-3444
Fax: (508) 540-8680
www.acciusa.com
Customer Service: custservice@acciusa.com
Technical Service: techservice@acciusa.com
Contract Test Service: testservice@acciusa.com

United Kingdom

Associates of Cape Cod Int'l., Inc.
Deacon Park, Moorgate Road
Knowsley, Merseyside, L33 7RX, United Kingdom
Tel: (44) 151-547-7444
Fax: (44) 151-547-7400
E-mail: info@acciuk.co.uk
www.acciuk.co.uk
UK Customer Service: customerservices@acciuk.co.uk
Company Registration Number: BR002906

Europe

Associates of Cape Cod Europe GmbH
Opelstrasse 14
D-64546 Mörfelden-Walldorf, Germany
Tel: (49) 61 05-96 10 0
Fax: (49) 61 05-96 10 15
E-mail: service@acciusa.de
www.acciusa.de

All Products and Services listed herein are offered exclusively under Associates of Cape Cod, Inc.'s Terms and Conditions of Sale which can be found at www.acciusa.com/pdfs/acctc.pdf



ASSOCIATES OF
CAPE COD
INCORPORATED

Specialists in Endotoxin and Glucan Detection

Corporate Headquarters

Associates of Cape Cod, Inc.

124 Bernard E. Saint Jean Drive
East Falmouth, MA 02536-4445 USA
Tel: (800) LAL-TEST (525-8378) or
(508) 540-3444
Fax: (508) 540-8680
www.acciusa.com
Customer Service:
custservice@acciusa.com
Technical Service:
techservice@acciusa.com
Contract Test Service:
testservice@acciusa.com

United Kingdom

Associates of Cape Cod Int'l., Inc.

Deacon Park, Moorgate Road
Knowsley, Liverpool, L33 7RX
United Kingdom
Tel: (44) 151-547-7444
Fax: (44) 151-547-7400
E-mail: info@acciuk.co.uk
www.acciuk.co.uk
UK Customer Support:
customerservices@acciuk.co.uk
Company Registration Number:
BR002906

Europe

Associates of Cape Cod Europe GmbH

Opelstrasse 14
D-64546 Mörfelden-Walldorf
Germany
Tel: (49) 61 05-96 10 0
Fax: (49) 61 05-96 10 15
E-mail: service@acciusa.de
www.acciusa.de